


Dr. Name _____ Phone # _____

Acct. # _____ Patient Name/ID _____
First Last

Address/Email _____ Deliver by 5 p.m. on _____



Tooth No. _____

Stump Shade _____

Final Shade _____

PONTIC DESIGN



ANTERIOR METAL DESIGN



POSTERIOR METAL DESIGN

- Full ceramic coverage
- Ceramic with lingual metal collar*
- 360° metal hairline or _____mm
- Metal occlusal excluding buccal cusp
- Metal occlusal including buccal cusp

INSTRUCTION FOR BUCCAL MARGIN

- Metal-ceramic junction margin*
- Ceramic butt margin

OCCUSAL STAINING

- None Light* Medium Dark

SCREW-RETAINED RESTORATIONS

- BruxZir Full-Strength* (w/ Ti-Base)
- IPS e.max (w/ Ti-Base)
- NEW!** BruxZir Esthetic (w/ Ti-Base)
- Obsidian Fused to White Noble
- Obsidian Fused to White High Noble

CERAMIC TO METAL

- Obsidian Fused to Non-Precious*
- Obsidian Fused to White Noble
- Obsidian Fused to White High Noble
- Yellow High Noble (76.6% Au)

ALL-CERAMIC/COMPOSITE

- BruxZir Full-Strength* (1,150 MPa)
- NEW!** BruxZir Esthetic (870 MPa) (stump shade recommended for restorations less than 1.5 mm thick)
- IPS e.max
- Bilayered CZ
- Generic Zirconia
- Composite
- Maryland Bridge

VIVANEER VENEER

- IPS e.max veneer*
- Layered IPS e.max veneer
- NEW!** BruxZir Esthetic veneer

PLAYSAFE MOUTHGUARDS

- Specify color(s) on Rx
- Upper* Lower
 - Light Pro Medium*
 - Heavy

FULL-CAST METAL

- Non-Precious
- White Noble (70.9% Pd)
- White High Noble (40% Au)
- Yellow Noble (41% Au)
- Yellow High Noble* (57.5% Au)
- Yellow High Noble (77% Au)

TEMPORARIES


- Hi-Tech Temps*
- Hi-Tech Temps with metal
- Diagnostic Wax-Up

BITE SPLINTS

- Upper* Lower
- Comfort H/S (hard, with soft reline)
- Soft Bite Splint

Enclosed with case: Impressions Models Bite Photos Other: _____ Try-In: Yes No

Rx SPECIFIC INSTRUCTIONS *NOTE: Please send a study model on all work involving anterior teeth.*



Implant System _____
 (if applicable)

Diameter _____
 (if applicable)

IF NO OCCLUSAL CLEARANCE

- Call doctor Spot opposing
- Metal occlusion Metal island
- Make this a permanent note in my master file

Signature _____

Submission of this Rx constitutes agreement with limited warranty terms and conditions. See reverse for details.

License # _____ Date _____

Place Shipping Airbill Tracking Sticker Here

CUSTOM ABUTMENTS

- Titanium Zirconia w/ Ti-Base
- Prepare existing abutment

TERMS AND WARRANTY INFORMATION

We honor VISA, MASTERCARD, AMEX and DISCOVER.

All rush cases must be prescheduled by calling **800-395-8205** before the case is shipped. Time of pickup and delivery may affect turnaround time.

TERMS: Cost of collection of any account will be paid by the customer. All accounts are payable within 30 days of statement date. **Accounts not paid within the stated terms will be subject to COD status and a late charge of 2 percent of the unpaid balance.** Prices subject to change without notice. Rx must be enclosed with original case submission.

NO-FAULT REMAKE POLICY: Smith-Sterling Dental Laboratories is pleased to process all remakes or adjustments at no additional charge if requested within the warranty period and accompanied by the return of the original appliance.



• BruxZir® Restorations



• Custom Abutments



• All-Ceramic Restorations
• PFM Restorations
• Full-Cast Restorations



• Nightguards
• Bite Splints
• Mouthguards



LIMITED WARRANTY/LIMITATION OF LIABILITY: For warranty terms and conditions and limitation of liability, [visit smithsterling.com/pura-vida-warranty](http://visit.smithsterling.com/pura-vida-warranty).




DENTISTS: There is no charge for one inbound and one outbound shipment per case. Additional shipments for die trims, bisque bake, coping or framework try-ins, reshades to a new shade, and oversized articulators will incur additional fees.

All restorations use *FDA-registered* dental alloys and ceramics. U.S. Food and Drug Administration Registration No. 3004540219


PREPARATION GUIDELINES

ANTERIOR: OBSIDIAN FUSED TO METAL



- A. 2.0 mm incisal reduction
- B. 1.5 mm middle third reduction
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

POSTERIOR: OBSIDIAN FUSED TO METAL



- A. 2.0 mm occlusal reduction
- B. 1.5 mm middle third reduction
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

CUSTOM ABUTMENTS AND SCREW-RETAINED CROWNS ARE COMPATIBLE WITH THE FOLLOWING IMPLANT SYSTEMS

BIOMET 3i™
Certain®

CAMLOG®
SCREW-LINE®

Dentium®
Implantium®†
SimpleLine® II†
SuperLine®†

DENTSPLY Implants
ANKYLOS® C/X*
ASTRA TECH Implant System®*
ASTRA TECH Implant System® EV*

Glidewell Direct
Hahn™ Tapered Implant System†
Inclusive® Tapered Implant System*

HIOSSEN®
HG System*

Keystone Dental
PrimaConnex®†

MegaGen
AnyRidge® Implant System*

Neoss®
Neoss® Implant System†

Nobel Biocare
Brånemark System® RP*
NobelActive®*
NobelReplace®*

Straumann®
Bone Level*
Tissue Level*

Sweden & Martina
Premium†
Shelta†

Zimmer Dental
Screw-Vent®*

*Manufactured using Inclusive® components. †Manufactured using original equipment manufacturer (OEM) components. Inclusive is a registered trademark of PrismaDent Dentalcraft, Inc. Hahn Tapered Implant is a trademark of PrismaDent Dentalcraft, Inc. All other trademarks are property of their respective owners. For Dentium implant systems, only titanium custom abutments are available.